SureFlex™ Laser Fiber

Single Use

Instructions for Use

SureFlex™ Laser Fiber <i>Single Use</i> Instructions for Use	page 1
Fibre laser SureFlex™ <i>à usage unique</i> Mode d'emploi	page 5
SureFlex™ Laserfaser <i>für den einmaligen Gebrauch</i> Gebrauchsanweisung	Seite 9
Fibra Laser SureFlex™ <i>Monouso</i> Istruzioni per l'uso	pagina 13
Fibra láser SureFlex™ <i>de aplicación única</i> Instrucciones de uso	página 17
Fibra óptica de laser SureFlex™ <i>Para uma única utilização</i> Instruções de utilização	pagina 21
SureFlex™ laserfiber <i>Eenmalig gebruik</i> Gebruiksaanwijsties	pagina 27

Read All Instructions Thoroughly Before Use.
Lire attentivement l'intégralité de ces instructions avant utilisation.
Alle Anweisungen vor dem Gebrauch sorgfältig durchlesen.
Leggere attentamente tutte le istruzioni prima dell'uso.
Lea atentamente todas las instrucciones antes de usarla.
Leia atentamente todas as instruções antes da utilização.
Lees alle instructies vóór gebruik aandachtig door.

Rx Only

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Instructions for Use

DESCRIPTION

The SureFlex™ laser fiber is a fiber optic laser delivery device consisting of a patented, high energy termination (Black Hole™), a length of silica/silica fiber with a UV cured polymer secondary cladding and an ethylene tetrafluoro-ethylene copolymer (ETFE) jacket. The standard fiber output is a defect-free laser polish tip. SureFlex fibers are designed for use in a wide variety of surgical procedures as an integral part of a surgical laser system.

INDICATIONS

SureFlex laser fibers are intended for use in laser-based surgical applications, including, but not limited to endoscopic, laparoscopic and open surgical procedures involving vaporization, ablation, coagulation, hemostasis, excision, resection and incision of soft and cartilaginous tissue and surgical procedures involving vaporization, ablation and fragmentation of calculi (urinary and biliary). SureFlex is designed primarily for holmium (Ho:YAG) lasers, but may be used with any laser wavelength between 500nm and 2200nm that has been cleared for surgical use, e.g. Nd:YAG, KTP, Alexandrite, Argon, Ruby, Diode. SureFlex laser fibers may be reprocessed within a single surgical session but may not be re-sterilized and reused.

INSTRUCTIONS FOR USE

1. Prior to activation of the laser, review the laser Operator's Manual for instructions for proper set-up and operation of the laser system.
2. Remove the fiber from the sterile pouch using aseptic technique. The SureFlex™ fiber is provided sterile (EO) and is not intended for reuse. Do not attempt to sterilize this product. Do not use the product if the pouch has been previously opened or if there is visual damage to the pouch.
3. Remove the fiber from the packaging starting from the distal end.
 → Do not start unwinding from the connector end.
 Unwinding from the connector end will make fiber removal from the packaging difficult and may damage the fiber.
4. Make sure the laser is "OFF" or in "STANDBY" mode before inserting the fiber into the laser aperture (SMA port). Remove the protective cap from the connector and inspect the glass ferrule at the laser launch face for contamination. It is advisable to wipe the input face clean with a low lint wipe prior to use to remove accumulated organic contaminants. Attach the connector to the laser and finger tighten the blue extension until a secure connection is achieved (generally indicated by release of lock-out mechanisms).

5. Handle the fiber with care as damage may occur if struck or bent sharply. Treat the output tip with particular care as it represents the most delicate, easily damaged portion of the assembly. If the fiber delivery system is damaged at the distal tip during handling or use, it may be repaired (below).
6. Activate the laser aiming beam and carefully inspect the full length of the fiber for breaks. Breaks will appear as bright spots along the length of the fiber. When the tip is held at about 5 cm from a light colored surface, the aiming beam should be visible as a uniform, round spot. If the spot appears dim it may be useful to slightly loosen the connector nut and rotate the shaft of the connector within the laser port to optimize the aiming beam appearance, particularly in the case of small core fibers. If no breaks are visible and the spot is uniform, the fiber is ready for use.
7. Beginning treatment: In “free beam” applications, hold the probe about 5 cm from the distal tip during use. If using an endoscopic delivery system, insert the probe until the tip extends approximately 1 cm from the distal end of the endoscope or is in the field of view. Adjust the laser power to the level appropriate for the intended use and fiber size:
6-watts maximum for SureFlex™ 150, 12-watts maximum for SureFlex™ 200, 35-watts maximum for SureFlex™ 273, 100 watts maximum for SureFlex™ 365, 550, & 910.
8. Intraoperative cleaning - If the tip accumulates debris during use, turn the laser to “STANDBY” then carefully wipe the tip with a sterile gauze pad. Do not scrape the fiber with sharp instruments, as damage may result. Debris that is permitted to accumulate and “char” on the fiber tip will accelerate degradation in performance and intraoperative reprocessing may be required.

INTRAOPERATIVE PROCESSING (DISTAL TIP ONLY)

INSPECTION

1. With a magnifying eyepiece (loupe), microscope, endoscope or video system (minimum magnification 20X), inspect the fiber face at the output to determine if processing is required. While inspecting the fiber, slightly tilt the fiber face to view a variety of angles and reflected light from the fiber surface.
2. The fiber faces should be free of any defects such as chips, cracks, scratches, pits, blemishes or debris within the central 90% of the output diameter (80% for SureFlex™ 150 and SureFlex™ 200). Some degradation of the fiber face at the working tip is normal and will occur with use. If pitting, cracks, chips or debris are observed, the fiber should be stripped and cleaved to restore the best performance and reduce the chances of a tip detachment.

CLEAVE AND STRIP THE DISTAL TIP OF THE FIBER

The recommended tools for cleaving and stripping fibers may be ordered from your authorized AMS distributor or directly from AMS. Check with an AMS sales representative for additional information. The following tools are recommended for stripping and cleaving the SureFlex™ laser fibers. See the accompanying

instructions packaged with the Fiber Cleaver and the Fiber Stripper for information on use (including sterilization).

FIBER SIZE	TOOL P/N	GUIDE	COLOR
150µm Autoclavable	80-10001-003	0.0135"	Lt. Blue
200µm Autoclavable	80-10002-003	0.018"	White
273µm Autoclavable	80-10003-003	0.018"	Dk. Green
365µm Autoclavable	80-10005-003	0.031"	Yellow
550µm Autoclavable	80-10011-003	0.046"	Black
910µm Autoclavable	80-10008-003	0.063"	Dk. Brown

Cleaving device for all fibers: 10-0760

NOTE: Verify the stripping tool corresponds with the fiber size (tools for other manufacturers' fibers and other sizes of SureFlex™ fibers look similar and are often confused). Use the blade and buffer guard retainer color as a guide. Insure tools are sterile and maintained within the sterile field for intraoperative fiber cleaving and stripping.

POTENTIAL COMPLICATIONS

Complications could include local and/or systemic infection, thermal damage to surrounding structures, local hematoma, dissection and perforation, tissue adhesion, distal tip detachment, and discomfort during and/or after energy application. In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of the tip. Properly reprocessed fibers (intraoperative), used within the accepted parameters of the field do not lose tips. If tips detach, review instructions.

CAUTIONS

- In soft tissues, increased laser exposure time will result in a deeper and wider zone of necrosis.
- Reprocessing and reuse of single use medical devices is hazardous and cannot be properly validated because the reprocessing facility cannot know what alterations may have been made to a product design after validation nor the extremes of materials properties that may not be adequately represented in a validation sample.
- This product should only be used by physicians who have been trained in proper use and operation of surgical laser systems and are familiar with the physiological applications of laser energy.
- If tissue appears to be stuck to the lasing tip, put the laser in standby and carefully remove the tip from the target tissue. Remove the fiber from the scope. Use a soft, sterile cloth to clean the debris from the glass tip, taking care

to avoid side pressure that may damage the tip.

- Do not bend fiber at sharp angles. If visible light (aiming beam) can be seen leaking from the fiber, fiber failure may result when therapeutic energy is applied as the fiber is deflected beyond the optical limits of total internal reflection.
- Care should be exercised with the glass tip to avoid severe impacts or side stresses that may fracture the tip.
- The time and energy required for treating a specific target area may be different when using the SureFlex™ laser fibers compared to using other fibers. A lower beam divergence, higher transmission efficiency and smaller laser spot size may alter treatment times.
- In soft tissue applications, during each lasing interval, it is helpful to move the distal tip back and forth or side to side 1 mm to 2 mm to prevent the tip from sticking to tissue.
- Immediately discontinue use if breaks or fractures appear in laser fiber. These breaks or fractures may allow undirected emission of laser energy, rendering the distal tip useless and potentially causing harm to surrounding tissues.
- **Always** insure all operating room personnel wear appropriate protective eyewear during the delivery of laser energy.
- Do **not** use the SureFlex™ laser fiber in the presence of flammable anesthetics or combustible materials.

HOW SUPPLIED AND STORAGE

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

- Sterility is guaranteed only if the package is not opened or damaged and the product is used within the expiration date on the package label.
- Store in a cool, dry place.

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SureFlex™ and Black Hole™ are trademarks of AMS.
All rights reserved. Printed in USA.
American Medical Systems Canada Inc.
U.S.A.
American Medical Systems Europe B.V.
© 2014 American Medical Systems, Inc.